

# ZEUS<sup>®</sup> Pipetting Module

## Specification Sheet



### Intelligent Liquid Handling


ZEUS (Z-Excursion Universal Sampler<sup>™</sup>) has been designed to take ownership of the entire pipetting process. It is based on Hamilton's revolutionary air displacement pipetting technology and uses disposable CO-RE II<sup>®</sup> tips to avoid carry-over.

With its integrated z-motion, ZEUS automatically picks up the tip, detects and follows the liquid level during aspiration and dispense. ZEUS adjusts its movements to the liquid type to ensure accurate and precise pipetting. The on-board Qualitative Pipetting Monitoring<sup>™</sup> (QPM) actively prevents malfunctions. These high level capabilities allow customers to significantly shorten development timelines and get product to market faster.


# ZEUS Product Drawings and Dimensions in mm

Straight Configuration


## Side




## Front



## Top



## Bottom



# Ordering Information

## ZEUS Pipetting Module

Part Number	Description	Includes
800700	ZEUS CO-RE module	
202977-0020	Test platform	Holder for the ZEUS module, and Communication box
202977-0021	Holder for ZEUS module	
202977-0022	Communication box	USB / Hamilton CAN converter, power supply, cables and software CD



Test platform

## Accessories and Spare Parts

Part Number	Description
800336	Holder for CO-RE tip rack
800337	Tip ejection sleeve (1/pk)
6606396-03	CO-RE II Tip Adapter (1/pk)



CO-RE II adapter



Tip ejection sleeve



Holder for CO-RE tip rack

## Disposable CO-RE II Tips\*

Part Number	Description				
235904-ZEUS	1000 µL tip	non-filtered	conductive	40 racks of 96	
235905-ZEUS	1000 µL tip	filtered	conductive	40 racks of 96	
235902-ZEUS	300 µL tip	non-filtered	conductive	60 racks of 96	
235903-ZEUS	300 µL tip	filtered	conductive	60 racks of 96	
235806-ZEUS	300 µL slim tip	non-filtered	conductive	40 racks of 96	
235647-ZEUS	300 µL slim tip	filtered	conductive	40 racks of 96	
235966-ZEUS	50 µL tip	non-filtered	conductive	60 racks of 96	
235948-ZEUS	50 µL tip	filtered	conductive	60 racks of 96	
235900-ZEUS	10 µL tip	non-filtered	conductive	60 racks of 96	
235901-ZEUS	10 µL tip	filtered	conductive	60 racks of 96	

\*Other tips are available on request



96 x 1000 µL CO-RE II tips

For more information on ZEUS and other Hamilton OEM Components please visit

[www.hamiltoncompany.com/OEM](http://www.hamiltoncompany.com/OEM)




## Product Specifications

### Dimensions


Width x Depth x Height

25 mm x 205 mm x 337 mm




Straight configurations  
(Single or parallel arrangement)

45 mm x 118 mm x 337 mm



Folded configurations  
(Single arrangement or as a pair)

161 mm x 102 mm x 337 mm



Right angle configuration  
(Single arrangement)

<b>Weight</b>	1100 g	
<b>z-Drive</b>	Maximum z-Excursion Resolution z-Axis Maximum Speed z-Axis	180 mm 0.1 mm 350 mm/s
<b>Pipetting Drive</b>	Volume Range Resolution Maximum Speed	1 µL – 1000 µL 0.1 µL 150 mm/s corresponding to 2'500 µL/s
<b>Pipetting Mode</b>	Single and aliquot dispensing	
<b>Process Security</b>	Tip presence sensor Pressure-based and capacitive Liquid Level Detection (cLLD and pLLD) Qualitative Pipetting Monitoring (QPM) Anti-Droplet Control (ADC)	
<b>Operating Data</b>	Power requirement Temperature range Relative humidity Noise level	48 VDC ± 5%, 1.8 A 18°C – 32°C 40% – 80% (non-condensing) < 60 dBA, Indoor use only
<b>Storage and Transportation</b>	Temperature range Relative humidity	-25°C to +70°C 40% – 80% (non-condensing, indoors)
<b>Communication</b>	CAN (Hamilton protocol), FCC connector Selectable baud rates Device addressing	125, 500 kBaud (others on request) up to 31 modules on a single CAN bus
<b>RoHS Compliant</b>	Yes	

## Pipetting Performance\*

Disposable Tip Size	Pipetting Volume	Accuracy	Precision
10 µL	1 µL	5.0%	4.0%
10 µL	5 µL	2.5%	1.5%
10 µL	10 µL	1.5%	1.0%
50 µL	1 µL	5.0%	4.0%
50 µL	5 µL	2.5%	1.5%
50 µL	50 µL	2.0%	0.75%
300 µL	10 µL	5.0%	2.0%
300 µL	50 µL	2.0%	0.75%
300 µL	300 µL	1.0%	0.75%
1000 µL	10 µL	7.5%	3.5%
1000 µL	100 µL	2.0%	0.75%
1000 µL	1000 µL	1.0%	0.75%

\*Pipetting specifications were determined gravimetrically using a high precision balance and strictly controlled environmental conditions: test temperature: 20 ± 2°C, relative humidity: 50% ± 5%, test liquid temperature range: ≤± 0.5°C of room temperature; balance used: Mettler Toledo MX5.

The measurements were done with Hamilton Verification Solution (deionized water with 0.1% NaCl and 0.01% Tween detergent) using standard CO-RE conductive tips. A new tip was used for each pipetting cycle (aspiration and dispense). At least 10 points were taken for each volume and pipetting module. Volumes > 20 µL were dispensed in jet mode. Volumes ≤ 20 µL were dispensed in surface mode. Results may vary using other liquid or environmental conditions.

© 2022 Hamilton Company. All rights reserved.  
All trademarks are owned and/or registered by Hamilton Company in the U.S. and/or other countries.  
Lit. No. 608770/06 — 07/2022

# HAMILTON®

Web: [www.hamiltoncompany.com](http://www.hamiltoncompany.com)

USA: 800-648-5950

Europe: +40-356-635-055

### Hamilton Americas & Pacific Rim

Hamilton Company Inc.  
4970 Energy Way  
Reno, Nevada 89502 USA  
Tel: +1-775-858-3000  
Fax: +1-775-856-7259  
[sales@hamiltoncompany.com](mailto:sales@hamiltoncompany.com)

### Hamilton Europe, Asia & Africa

Hamilton Central Europe S.R.L.  
str. Hamilton no. 2-4  
307210 Giarmata, Romania  
Tel: +40-356-635-055  
[contact.lab.ro@hamilton-ce.com](mailto:contact.lab.ro@hamilton-ce.com)

To find a representative in your area, please visit [hamiltoncompany.com](http://hamiltoncompany.com).